

**About JINNI MD<sup>®</sup>:** JINNI MD<sup>®</sup> will help promote blood circulation, ease discomfort of fatigue, swelling & varicose veins.

**JINNI MD<sup>®</sup> Applications & Uses:** Relief from tired and aching veins, helps to lessen appearance of varicose veins, reduce mild swelling, lower risk of DVT, anti-embolism, anti-thrombosis & makes your legs feel healthy.

JINNI MD<sup>®</sup> is available in a full range of sizes, compressions, colours and Open-Toe or Closed-Toe styles to suit all individuals.

### What size is right for me?


\*To determine what size JINNI MD<sup>®</sup> COMPRESSION STOCKING is required, it is important to measure the following:

- (cB) Smallest ankle circumference
- (cC) Largest calf circumference
- (cD) Circumference below the knee
- (cF) Largest thigh circumference
- (cG) Upper thigh circumference

\*All measurements should be completed first thing in the morning before active movement is undertaken.

\*Each leg should be measured individually & the lengths should be measured (ID, IF, IG) as required.

*REFER TO Following SIZE CHART*



### What compression should I get?

#### Therapeutic support (15-20 mmHg Mild Compression)

Help relief from the tired, aching, fatigued legs, minor ankle and leg swelling, minor varicosities, incipient varicosis during pregnancy, prevention of thrombosis.

#### Class I (20-30 mmHg Moderate Compression)

Help relieve moderate ankle and leg swelling, moderate varicosities, post-sclerotherapy applications.

#### Class II (30-40 mmHg High Compression)

For chronic venous insufficiency, pronounced oedema, atrophie blanche, dermatosclerosis and after resolution of severe ulcers.

### What style?

- Knee High (Open/Closed Toe)
- Thigh High (Open/Closed Toe)
- Waist High (Open/Closed Toe)
- Sporting/Travel Socks

TYPE	KNEE HIGH STOCKINGS			
	Circumference measurements (cm)			Length (cm)
S	19 - 21	29 - 34	25 - 31	38 - 39
M	21 - 23	33 - 37	30 - 36	39 - 40
L	23 - 25	36 - 41	35 - 40	40 - 41
XL	25 - 27	40 - 45	38 - 45	41 - 42
XXL	27 - 29	44 - 47	44 - 48	42 - 43

TYPE	THIGH HIGH/PANTYHOSE STOCKINGS			
	Circumference measurements (cm)		Length (cm)	
S	38 - 45	47 - 55	58 - 59	68 - 71
M	44 - 49	54 - 62	59 - 60	71 - 72
L	48 - 53	61 - 66	60 - 61	72 - 73
XL	52 - 57	65 - 70	61 - 64	73 - 75
XXL	56 - 61	69 - 75	64 - 66	75 - 78

Thigh/waist high stockings: Measure cB,cC,cD,cF,cG, and the lengths.

### Fitting Recommendations

Before wearing compression stockings, it is important that the correct size is determined.

- Narrowest ankle circumference
- Check against measuring guide
- Match measurements with the required amount of compression

It is important to be shown EARLY the correct way to fit & wear compression stockings/socks.

Comfort & effectiveness will increase the patients' willingness to wear the stockings.

When putting on stockings, patients should remove all jewellery and try to avoid snagging the material with their fingernails.

If size/compression/style is out of stock, see shop staff to have it ordered in.

#### Disclaimer:

Always read the instructions/label.  
Use only as directed.

#### Global Distributor

WallCann Pty Ltd  
147 – 149 Holbrooks Rd  
Underdale SA 5032  
Ph: +61 8 8352 1006  
www.jinnimd.com

**TGA**  
APPROVED